

DE LA RECHERCHE À L'INDUSTRIE



PathoTRACK platform on IFB Cloud

IFB Cyclone Hackathon

Adrien JOSSO

IFB-Core - Genoscope/LABGeM

CBRNE Research and Innovation Project

- ▶ Partners
 - ▶ CEA - CNS / LABGeM
 - ▶ CEA - CNG / Sequencing Platform
 - ▶ IFB-Core
 - ▶ DGA - Maîtrise NRBC
- ▶ started since 2015
- ▶ born from CBRNE R&I project named *SynBioWatch* (2012-2014)



Wet Lab / Sequencing


- ▶ Improve DNA extraction protocol (complex matrix: soil, petrol, ...)
- ▶ Benchmarck "**portable**" sequencing technologies
 - ▶ Long reads: Oxford Nanopore Minion, IonTorrent PGM
 - ▶ *versus*
 - ▶ Short reads: Illumina MiSeq

Bioinformatics

- ▶ **Quick identification** and **characterization** of pathogens from environmental and clinical samples
- ▶ Develop reference databases (Taxonomique Assignment, ...)
- ▶ Bioinformatic tools integration (SortMeRna, Diamond, Kraken, Centrifuge)
- ▶ Workflow management system: SnakeMake


PathoTRACK platform

Workflow Management System



PathoTRACK platform

Analysis




3 workflows

- ▶ Preprocessing sequencing data
- ▶ Functional annotation (genes involved in drug resistance and virulence)
- ▶ Taxonomic assignment

PathoTRACK platform

On IFB Cloud



Micro-Services Architecture

- ▶ 4 services
 - ▶ Web User Interface
 - ▶ LIMS
 - ▶ Workflow Manager
 - ▶ MySQL server
 - ▶ Each service is elastic, resilient, composable, minimal, and complete
-
- ▶ OS: Linux CentOS 6 (7 in progress) (64 bits)
 - ▶ Can be deployed on Desktop, VM, servers, HPC clusters, Cloud (IFB)
 - ▶ Start and stop the platform using a daemon

Installation

- ▶ The platform is packaged in an archive (tar.gz)
- ▶ Bash script (need root privileges)
- ▶ Update the system
- ▶ Install all needed packages using yum (they will be packaged with the platform to perform an offline installation)
- ▶ Install all services of PathoTRACK platform
- ▶ Configure the automatic launch of the platform (init.d / chkconfig)
- ▶ Configure Apache httpd server in ProxyPass mode (ports redirecting)
- ▶ Initialize the platform

Use the CYCLONE components to deploy the platform on the new version of IFB cloud

- ▶ Use SlipStream to deploy the platform ? each service on a dedicated VM ?
- ▶ Configure VM access (Authentication)
- ▶ Configure a dedicated cluster on IFB Cloud ?



Thank for your attention

LABGeM

Project Leader: Stephane Cruveiller, Ph D.

- ▶ Johann Beghain
- ▶ Marine Séjourné

CNG Sequencing platform

- ▶ Robert Olasso, Ph D.
- ▶ Céline Bigot, Ph D.

IFB-Core

- ▶ Adrien Josso